BLOODLESS PLUS CAS10

UNA AUTENTICA REVOLUCION EN LA CIRUGIA
A REAL REVOLUTION IN THE SURGICAL FIELD

INCASMEDICAL

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

Advantages and

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Foto 3 Y 4 dano termico difundido

Removal of Capillar Hemangioma and Lobular Hemangioma with wire loop Bloodless Plus

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Clinical Anatomical study on the employment in phatology and Dermosurgery and histological-clinical study of the pulmonary tissues using Bloodless Plus

Bloodless Plus is a Molecular Resonance Lancet using a new physical principle in order to carry out cutting and coagulation at under 50°C & 65°C temperature respectively within biological tissues. This new physical principle exploits properly modulated high frequency waves in order to obtain the separation of the tissue acting on the molecular ties of the cells and not by thermal vaporization. The obtained result is therefore a cold cut (under 50°C) and the tissue to the edges of the cut are without eschar, with the advantage to reduce the danger of infection, post-operation pains and recovery time.

It generates the frequencies 4/8/12 MHz or more (any power/function) within cell safety spectrum derived from molecular resonance principle that too in resonance with the frequencies of molecular cellular bonds.

Rapid biopsies can be performed, because the samples taken are not damaged, with consequent reliable analysis.

To prove what has just been written above, there are some samples taken in dermatology by **Doctor Massimo Cogo**, of the Clinic "Villa Berica" of Vicenza and analyzed by Doctor Raffaele Giannotti, specialist in pathologic anatomy and histopathology at the 1st Dermatological clinic of the University in Milan.

We also show some examples of pulmonary tissues taken by **Prof. Federico Rea** at the Thoracic Surgery Department at the University of Padova (Responsible) **Prof. Francesco Sartori**, and analysed by Doctor Fiorella Calabrese in the histological laboratory (Responsible: **Prof. Marialuisa Valente**).

The fragment of pulmonary parenchyma is cut by a molecular resonance lancet Bloodless Plus.

As you can see, there are no alterations around the cut neither in the histological (A) nor in the ultrastructural (B) preparations.

The fragment of pulmonary parenchyma (pictures C and D) lacking the pleural area proves to be treated by the molecular resonance lancet Bloodless Plus, in order to obtain a clot. As you can see, neither in the histological (C) nor in the ultrastructural (D) preparations, the alveolus is damaged.

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Therapeutic applications : Bloodless Plus only needs a special module for its ready use in therapeutic applications like Edema (post-operatively or post-traumatic), Neck Pain, Back Pain, Muscle Injury, Sprain, Arthritis, Osteoarthritis, Disc Herniation and many more. The hand-piece with electrodes up to 5 cm dia can be connected or patches or adhesive electrodes can be applied with the help of user friendly programs. During therapeutic applications mode the maximum output of 75 watt at 50Ω to 100 watt at 1000Ω can be derived.

The frequencies of cell safety spectrum are transmitted to the patient's area of concern to induce anti-inflammatory and anti edema effect resulting into substantial decrease in pain. With the continuous session the damaged tissue regenerates in a stable form without any cellular damage and side effects. Bloodless Plus has been proved a useful tool in sports medicine also.

Bloodless Plus: a universal instrument


Bloodless Plus is an instrument suitable for different surgical fields like urology, gynaecology and orthopaedics.

In particular, it was created for a use in ENT, General Surgery, Neuro Surgery, Plastic Surgery, Dermatology and many more. Bloodless Plus is extremely easy to use and its users don't need special trainings.

A very important feature of Bloodless Plus is constituted by the possibility of creating every kind of tool for cut and clot following the instructions of the surgeon that needs it for special operations including endoscopic approach for the treatment of Degenerative Disc Diseases.

It offers a wide range of accessories; lancets of different sizes, spatula, loops with a lot of shapes, bipolar forceps for blood vessels' clot and tools for clot, which make of it a really complete and essential instrument in the surgical field.

Moreover, Bloodless Plus well matches with all kind of equipment of the operating room, defibrillator included, thanks to the technology by which it has been manufactured.

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We request our channel partners and customers to get in touch with us by email only

Certification & Inspection
ISO 9001:2015 Cert. No. UQ-103041816
ISO 13485:2016 Cert. No. UQ-103041817